




ATTENTION – BEFORE APPLICATION:

- 1) **PREPARE:** Make sure you have Ivory or Gold Dial bar soap – no liquid soap. Consider using a Kind Klamp™ if needed (see below).
- 2) **REVIEW:** Familiarize yourself with the product and review the application instructions ([see back](#)).
- 3) **WATCH:** Watch the application video at www.MensLiberty.com.
- 4) **CALL:** Call our Nurse Support Line for an appointment **BEFORE** applying Men's Liberty for the first time.

**Please allow up to 24 hours for a return call.*



NURSE SUPPORT LINE

(888) 920-9363

Hours: M-F, 9:00 AM-5:00 PM Eastern

Manufactured by:

BIODERM®

COMPASSIONATE SOLUTIONS

Removing Men's Liberty

FreeDerm™ is recommended to remove Men's Liberty. Men's Liberty can also be removed by soaking the seal in warm water or a wet washcloth until it is milky white and rolls gently from the skin.

TO USE: Spray the edge of the adhesive, lifting gently from the skin as you go. Continue to spray under the adhesive until completely removed.

Never pull or tear from the skin.



The Kind Klamp Penile Clamp

Kind Klamp is designed to temporarily stop the flow of urine. Consider ordering a Kind Klamp if you think you may experience leaking during application, or if you are uncircumcised or have retraction.

Call us at (800) 373-7006 to order.



BioDerm, Inc.
 12320 73rd Court North
 Largo, FL, 33773, USA
 Orders: 800-373-7006
 Support: 800-864-6755
 CustomerCare@BioDermInc.com
 www.BioDermInc.com



Authorized Representative:
 Advna Ltd.
 Pure Offices, Plato Close,
 Warwick, CV34 6WE, UK

**NOT MADE WITH
 NATURAL RUBBER LATEX**

©2017 BioDerm, Inc. All rights reserved.
 BioDerm® is a registered trademark of BioDerm, Inc. Men's Liberty™, FreeDerm™, BioPlus+™ and Kind Klamp™ are trademarks of BioDerm, Inc.
 Men's Liberty US patent #7,896,857, #8,551,062.

1 Begin by washing the head of the penis with Ivory or Gold Bar Dial soap. Allow skin to dry completely. Use a Kind Klamp if needed.



2 Wipe the head of the penis with the BioPlus+™ wipe. Allow 15-20 seconds to dry until the skin is no longer tacky.



! About Your First Few Days


As a result of using pads or condom catheters, your skin has **A LOT** of moisture. Men's Liberty **LOVES** moisture and will soak it all up until your skin returns to a normal moisture balance. It can take up to three consecutive applications for you to get a full 24-hour wear time. **Always reapply another Men's Liberty right away** – don't go back to other products as they will reintroduce moisture!

3 Position Men's Liberty so that the vent is facing up. Remove release papers and attach the faceplate to the tip of the penis with the center opening directly in line with the urinary opening. Smooth down petals.



4 Apply the faceplate seal around the head of the penis, flat and wrinkle free. Be sure the hydrocolloid is pressed into all the ridges and grooves in the skin. Press and hold for 10-15 seconds to secure the seal.

Flip over for removal instructions.



For overnight use: Unplug the bottom of the Men's Liberty pouch and connect your bedside bag. Hang the bed bag from your bed frame or rest against the floor. Keep below the bladder to ensure drainage. Do not twist or kink the tube.

If irritation occurs, discontinue use and consult your healthcare professional.