

Men's Liberty™ Instructions for Use

We are here to support you!
Call the Nurse Support Line : (888) 920-9363
Hours: Mon-Fri, 9:00am-5:00pm EST

Follow these steps for proper application:

1 Review Contents

- Familiarize yourself with the products to the right and review the application instructions (see reverse side) before applying.
- Make sure you have a non-moisturizing bar soap like Ivory or Gold Dial at home - do not use liquid soap.



- Consider using a Kind Klamp™ if you think you may experience leaking during application, you are uncircumcised or have retraction.

Call (800) 373-7006 to place an order.

- Watch the application instruction video at mensliberty.com

Don't Run Out!

Call your distributor to place a re-order when you have 10 units or less.

Preparation



BioPlus™ Prep Wipe

Optional Leak Prevention



Kind Klamp™
Penile Clamp

Removal



FreeDerm™
Adhesive Remover

Application



Faceplate

Circular Tube

Vent

Men's Liberty™
Faceplate Strip

Urine Collection
Pouch

Drain Plug



Bioderm, Inc.
12320 73rd Court North
Largo, FL 33773, USA
Orders: (800) 373-7006
Support: (800) 920-9363
CustomerCare@MensLiberty.com
www.MensLiberty.com

EC REP CE
Authorized Representative:
Advena Ltd.
Pure Offices, Plato Close,
Warwick, CV34 6WE, UK

**NOT MADE
WITH NATURAL
RUBBER LATEX**

©2018 BioDerm, Inc. All rights reserved. BioDerm® is a registered trademark of BioDerm, Inc. Men's Liberty™, FreeDerm™, BioPlus+™ and KindKlamp™ are trademarks of BioDerm, Inc. Men's Liberty US patent #8,551,062, ML US 7,896,857, Penile Clamp US Patent: 767761. Additional US and International patents are pending.

PART# 23046 LT 95207 A

Men's Liberty™ Instructions for Use

2 Wash & Dry



Begin by washing the head of the penis with a non-moisturizing bar soap like Ivory® or Gold Dial®. Dry thoroughly with a cotton or paper towel. Use a Kind Klamp to prevent leaking during washing if necessary. Dry thoroughly.

3 Pre-Application Prep



Wipe the head of the penis with BioPlus+™ wipe. Allow 15-20 seconds to dry until the skin is no longer tacky.

4 Positioning & Application



Position Men's Liberty so that the vent is facing up. Align the center opening of the circular tube and petals with your urinary opening (meatus). Remove release papers around the petals and attach the faceplate to the tip of the penis with the center opening directly in line with the urinary opening. Smooth down petals.

5 Secure and Seal



Position the faceplate strip under the tip of the penis and center in the middle. Remove the adhesive backing and ensure the edge is over both the end of the clear circular tube and the tip of the penis. Wrap the strip around each side, making sure it is flat and wrinkle free. Use your fingers to push the ridges and grooves so there are no air pockets. Press and hold for 10-15 seconds to secure the seal.

6 Removal



FreeDerm™ is recommended to remove Men's Liberty. Men's Liberty can also be removed by soaking the seal with warm water or a wet washcloth until it is milky white and rolls gently from the skin. **TO USE:** Spray the edge of the adhesive, lifting gently from the skin as you go. Continue to spray under the adhesive until completely removed. **Never pull or tear from the skin.**







Your First Few Days

As a result of using diapers or condom catheters, your skin may have a lot of moisture. Men's Liberty will soak up extra moisture until your skin returns to a normal balance. It can take up to three consecutive applications for you to get a full 24-Hour wear time. Always reapply another Men's Liberty right away if one falls off. It may take up to a full week of wear-time until you are comfortable wearing Men's Liberty.

Nurse Support Line: **(888) 920-9363**

Hours: Mon-Fri, 9:00am-5:00pm EST



For Overnight Use

Unplug the bottom of the Men's Liberty pouch and connect your bedside bag. Hang the bed bag from your bed frame. Keep below the bladder or bed to ensure drainage. Do not twist or kink the tube.

If you have any questions, call the nurse support line and we will walk you through each step of your first application. **(888) 920-9363**